



[SHOP NOW](#)

How to take care of your Ribcap?

Here are some useful tips for maintaining your Ribcap. If you have any questions or concerns, feel free to reach out to our dedicated [support team](#). We're here to assist you in keeping your Ribcaps in excellent condition!



Hand Washing Method:

1. ***Prepare the Basin:***
 - Start by filling a basin or sink with lukewarm water.
 - Add a mild detergent, steering clear of harsh chemicals like bleach.
 2. ***Submerge and Agitate:***
 - Gently submerge your Ribcaps into the soapy water.
 - Give it a gentle agitating motion to help remove dirt and stains effectively.
 3. ***Soak and Rinse:***
 - Allow your Ribcaps to soak for a few minutes, letting the detergent work its magic.
 - After soaking, carefully rinse your Ribcaps under cool, running water to remove all detergent residues.
-

Spot Cleaning Method:



1. *Soapy Sponge or Cloth:*

- Opt for the spot cleaning method by using a sponge or cloth soaked in soapy water.

2. *Gentle Scrubbing:*

- Gently scrub the soiled areas with the soapy sponge or cloth, targeting specific stains.

3. *Rinse with Care:*

- Rinse the cleaned areas with clean water, ensuring no soap residues remain.
-



Drying:

1. *Lay Flat and Pat Dry:*

- Lay your Ribcaps flat on a clean towel.

- Gently pat them to remove excess water, being mindful of the delicate materials.

2. *Air Dry Naturally:*

- Allow your Ribcaps to air dry naturally. Avoid exposing them to direct sunlight or artificial heat sources.

Important Reminder: Do Not Machine Wash:

- It's crucial to avoid machine washing or tumble drying your Ribcaps, as this could potentially damage the materials.

Can I use a Ribcap with hearing aids?

Ribcap's design focuses on providing protection for the head, particularly against impacts and abrasions. While wearing a Ribcap, you should be able to use a hearing aid simultaneously. The material is typically flexible and should not interfere with the placement or function of a hearing aid.



See what our customers think about our products.



I get many compliments on my helmet and i highly recommend ribcap to anyone who is in need of a helmet.

Joann R.



Recommend highly as the head protection i wear in the house is not something i would want to wear in public.

Roxanne T.



The hat is warm, but my ears get cold so i can't wear in really cold weather.

Whether you're suffering from a disorder like epilepsy, autism, cerebral palsy, or Parkinson's disease, or simply want to enhance your quality of life, Ribcap was designed for everyone - including active children, young adults, and seniors facing all sorts of life challenges!

No longer want to receive these emails? [Unsubscribe.](#)

Ribcap Zwaarveld 34 | 9220 Hamme | Belgium Hamme, Hamme 9220

